

Name : Vaccine Carrier
Brand : BISAMED
Type : BSM-VC15



Features :
- Capacity 15 L
- Electronic Control
- Foam Material Polyurethane Cycloisopentane
- Phase Change Material (PCM) Cold Pack
- Battery for Electronic Control and Emergency Cooling

Specifications :
Temperature : 2 - 6 °C
Capacity : 15L
Power Supply : AC 110 - 240V AC/ 12-24V DC
Material : Foam Polyurethane Cycloisopentane
Battery : For Electronic Control and Emergency Cooling
Cold Protection : Passive Cooling Mode Advance Design