

Name : Vaccine Carrier
Brand : BISAMED
Type : BSM-VC50



Features :

- Capacity 50 L
- Electronic Control
- Foam Material Polyurethane Cycloisopentane
- Phase Change Material (PCM) Cold Pack
- Battery for Electronic Control and Emergency Cooling

Specifications :

Temperature : 2 - 6 °C
Capacity : 50L
Power Supply : AC 110 - 240V AC/ 12-24V DC
Material : Foam Polyurethane Cycloisopentane
Battery : For Electronic Control and Emergency Cooling
Cold Protection : Cooling System Compressor DC