

Name : Extraoral Aerosol Dental Suction
Brand : BISAMED
Type : BSM-ADS01



Features :
- Low Noise to make operator more comfortable
- Touchless button On/Off avoiding cross contamination.
- High speed and pressure motor pump to increase the performance

Specifications :

Hepa Filter (H13) : 0,3 mikron with Carbon & Nano AG+
UV-C : 6 Watt, 254 nm Wavelength
Pump Motor : AC 220V / 50 H
Power Consumption : Max. 572 watt power consumption
Speed : Approx. 0 - 22800 rpm
Pressure : Approx. 11.000 Pa
Noise : Approx. 55 - 85 dB
Button On/Off : Touchless