



Exclusive Distributor :
PT. BINTANG MONO INDONESIA



PROFILE

PT. Bintang MONO Indonesia adalah importir dan distributor utama produk alat kesehatan In Vitro Diagnostic dan Cardiovascular Medical Device bersertifikasi CDAKB KEMENKES RI.

Kami menyediakan produk berkualitas dan membangun reputasi sebagai distributor terkemuka dibidang alat kesehatan.

Kami memiliki pengalaman panjang dan terpercaya dalam industri alat kesehatan dengan menyediakan alat kesehatan yang aman dan berkualitas tinggi untuk kami tawarkan kepada sub-distributor dan konsumen dengan pilihan harga terbaik.

Contact

PHONE:

+62 21 29040042

WEBSITE:

www.bintangmonoindonesia.com

PRODUK-PRODUK PENANGGULANGAN PANDEMI CORONAVIRUS



Coronavirus **ANTIGEN** Rapid Test Cassette (Swab)



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand Internasional
- ☒ Distributor bersertifikasi CDAKB

Link barcode produk:



Product Information

SPESIFICATION

Product Name

Healgen coronavirus Ag Rapid Test Cassette (swab)

Reg. Number

KEMENKES RI AKL 20303026076

Specimen

Naso-Nasal Swab

Package

20 Test/Kit

Type

Cassette

Function

SARS viral antigen is generally detectable in upper respiratory specimens during the acute phase of infection. The Coronavirus Ag Rapid Test is a rapid screening tool to detect the presence of SARS viral antigen in the form of a visually interpreted result within minutes.

PRODUK-PRODUK PENANGGULANGAN PANDEMI CORONAVIRUS



COVID-19 IgG/IgM Rapid Test (whole blood/serum/plasma)



CFDA

FDA



TGA



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand Internasional
- ☒ Distributor bersertifikasi CDAKB

Link barcode produk:



SPESIFICATION

Product Name

HEALGEN COVID-19 IgG/IgM Rapid Test Cassette (Whole Blood/serum/plasma)

Reg. Number

KEMENKES RI AKL 20303023810

Specimen

Whole blood/serum/plasma)

Package

25 Test/Kit

Type

Cassette

Function

COVID-19 IgG/IgM Rapid Test Cassette(Whole Blood/Serum/Plasma) is a solid phase immunochromatographic assay for the rapid, qualitative and differential detection of IgG and IgM antibodies to 2019 Novel Coronavirus in human whole blood, serum or plasma.

**PRODUK-PRODUK PENANGGULANGAN
PANDEMI CORONAVIRUS**



**VTM Non-Inactivated
Double Swab**



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand Internasional
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



Product Information

SPESIFICATION

Product Name

RONGYE Viral Transport Tube with Specimen Colloection Swab (Non-Inactivated Preservative Medium)

Reg. Number

KEMENKES RI AKL 10302023343

Specimen

Oral Swab/Nasopharyngeal Swab

Package

40 Test/Kit

Type

Cassette

Function

The viral preservative medium could be used for collecting, preserving and transporting the specimens of flu, Novel Coronavirus , Bird Flu, Hand-foot-mouth Virus, Morbilli cirus and other nasopharynx pathogen

**PRODUK-PRODUK PENANGGULANGAN
PANDEMI CORONAVIRUS**



**VTM Inactivated
Double Swab**



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand Internasional
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPESIFICATION

Product Name

RONGYE Viral Transport Tube with Specimen Colloection Swab (Inactivated Preservative Medium)

Reg. Number

KEMENKES RI AKL 10302022428

Specimen

Oral Swab/Nasopharyngeal Swab

Package

40 Test/Kit

Type

Cassette

Function

The viral preservative medium could be used for collecting, preserving and transporting the specimens of flu, Novel Coronavirus , Bird Flu, Hand-foot-mouth Virus, Morbilli cirus and other nasopharynx pathogen

RAPID TEST DRUG OF ABUSE



Drug Of Abuse Multi 6 in 1

(AMP, THC, MOP, COC, MET, BZO)



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand Internasional
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPEISIFICATION

Product Name

HEALGEN DOA Multi Panel 6 in 1
(AMP, THC, MOP, COC, MET, BZO)

Reg. Number

KEMENKES RI AKL 20103610974

Specimen

Urine

Reading Time

5 - 10 Menit

Package

50 Test/Kit

Type

Panel

Parameter

Amphetamines (AMP), Morphine (MOP/MOR)
Marijuana (THC), Cocaine (COC)
Methamphetamine (MET), Benzodiazepine (BZO)

Function

Detect the presence of 6 types of drugs in urine containing narcotics as well as psychotropic.

RAPID TEST DRUG OF ABUSE



Drug Of Abuse Multi 5 in 1

(AMP, THC, MOP, COC, MET)



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand Internasional
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPEISIFICATION

Product Name

HEALGEN Multi Drug og Abuse Multi 5
(AMP, THC, MOP, COC, MET)

Reg. Number

KEMENKES RI AKL 20103718058

Specimen

Urine

Reading Time

5 - 10 Menit

Package

50 Test/Kit

Type

Panel

Parameter

Amphetamines (AMP), Morphine (MOP/MOR),
Marijuana (THC), Cocaine (COC),
Methamphetamine (MET)

Function

Detect the presence of 5 types of drugs in urine containing narcotics as well as psychotropic.

RAPID TEST DRUG OF ABUSE

HEALGEN
Drug Of Abuse Multi 3 in 1
(MOP, AMP, THC)

e catalogue



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand International
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPESIFICATION

Product Name

HEALGEN Multi Drug of Abuse Multi 3 (MOP, AMP, THC)

Reg. Number

KEMENKES RI AKL 20103614262

Specimen

Urine

Reading Time

5 - 10 Menit

Package

50 Test/Kit

Type

Panel

Parameter

Morphine (MOP), Amphetamines (AMP), Marijuana (THC)

Function

Detect the presence of 3 types of drugs in urine containing narcotics as well as psychotropic.

RAPID TEST INFECTIOUS DISEASE

HEALGEN
HCV Test Cassette

e catalogue



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand International
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPESIFICATION

Product Name

HEALGEN HCV Test Cassette

Reg. Number

KEMENKES RI AKL 30305919146

Specimen

Serum / Plasma

Reading Time

15 Menit

Package

50 Test/Kit

Type

Cassette

Sensitivity

98,13 %

Spesificity

98,9 %

Function

The qualitative detection of antibody to Hepatitis C Virus in serum or plasma

RAPID TEST INFECTIOUS DISEASE



HBsAb Test Cassette



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand International
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPESIFICATION

Product Name

HEALGEN HBsAb One Step Rapid Test Cassette

Reg. Number

KEMENKES RI AKL 30305020010

Specimen

Serum / Plasma

Reading Time

15 Menit

Package

50 Test/Kit

Type

Cassette

Sensitivity

99,0 %

Spesificity

98,5 %

Function

The qualitative detection of Hepatitis B Surface Antibody in serum or plasma.

RAPID TEST INFECTIOUS DISEASE



HBsAg Test Cassette



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand International
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPESIFICATION

Product Name

HEALGEN HBsAg Test Cassette

Reg. Number

KEMENKES RI AKL 30305919153

Specimen

Serum / Plasma

Reading Time

15 Menit

Package

50 Test/Kit

Type

Cassette

Sensitivity

99,4 %

Spesificity

99,5 %

Function

The qualitative detection of Hepatitis B Surface Antigen in serum or plasma.

RAPID TEST INFECTIOUS DISEASE



Syphilis Test Cassette



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand Internasional
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



Product Information

SPESIFIKATION

Product Name

HEALGEN Syphilis Test Cassette

Reg. Number

KEMENKES RI AKL 20303612283

Specimen

Serum / Plasma

Reading Time

15 Menit

Package

50 Test/Kit

Type

Cassette

Sensitivity

98,1 %

Spesificity

98,9 %

Function

The qualitative detection of antibodies (IgG and IgM) to Treponema Pallidum (Syphilis) in serum or plasma to aid in the diagnosis of Syphilis.

RAPID TEST TROPICAL DISEASE



Typhoid Test Cassette



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand Internasional
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPESIFIKATION

Product Name

HEALGEN Typhoid Test Cassette

Reg. Number

KEMENKES RI AKL 20303717667

Specimen

Serum / Plasma

Reading Time

10 - 15 Menit

Package

50 Test/Kit

Type

Cassette

Sensitivity

IgM = 91,8 % , IgG = 92,7 %

Spesificity

IgM = 98,7 % , IgG = 98,8 %

Function

The simultaneous detection and differentiation of anti-Salmonella typhi (S. typhi) IgG and IgM in human serum and plasma.

RAPID TEST TROPICAL DISEASE

HEALGEN Malaria (Pf/Pv) Test Cassette

e catalogue

- ✓ Terdaftar di E-CATALOGUE LKPP
- ✓ Teregistrasi KEMENKES RI AKL
- ✓ Bersertifikasi Internasional
- ✓ Brand Internasional
- ✓ Distributor bersertifikasi CDAKB



Link barcode produk:



SPESIFICATION

Product Name

HEALGEN Malaria (Pf/Pv) Test Cassette

Reg. Number

KEMENKES RI AKL 20303613230

Specimen

Whole Blood

Reading Time

20 - 30 Menit

Package

50 Test/Kit

Type

Cassette

Sensitivity

(Pf+) = 100 % (Pv+) = 100 %

Specificity

(Pv+) = 98,7 % (Pv+) = 98,7 %

Function

The Cassette simultaneous detection and differentiation of Plasmodium falciparum (Pf) and vivax (Pv) antigen in human blood specimen.

RAPID TEST TROPICAL DISEASE

HEALGEN Dengue (IgG/IgM) Rapid Test Cassette

e catalogue

- ✓ Terdaftar di E-CATALOGUE LKPP
- ✓ Teregistrasi KEMENKES RI AKL
- ✓ Bersertifikasi Internasional
- ✓ Brand Internasional
- ✓ Distributor bersertifikasi CDAKB



Link barcode produk:



SPESIFICATION

Product Name

HEALGEN Dengue IgG/IgM One Step Rapid Test Cassette

Reg. Number

KEMENKES RI AKL 20303615192

Specimen

Whole Blood/Serum/Plasma

Reading Time

10 - 15 Menit

Package

50 Test/Kit

Type

Cassette

Sensitivity

95,8 %

Specificity

99,0 %

Function

Test (Cassette) The qualitative detection of IgG and IgM antibodies to Dengue virus in human whole blood, serum, or plasma as an aid in the diagnosis of primary and secondary Dengue infections.

RAPID TEST INFECTIOUS DISEASE



TB Test Cassette



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand Internasional
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPEISIFICATION

Product Name

HEALGEN TB Test Cassette

Reg. Number

KEMENKES RI AKL 20303021570

Specimen

Serum/Plasma

Reading Time

15 Menit

Package

50 Test/Kit

Type

Cassette

Sensitivity

IgM = 83,5 % , IgG = 87,2 %

Spesificity

IgM = 98,1 % , IgG = 97,3 %

Function

The rapid detection of IgM anti-Myco-bacterium Tuberculosis and IgG anti-Myco-bacterium Tuberculosis in human serum or plasma.

RAPID TEST INFECTIOUS DISEASE



HIV 1/2 Test Cassette



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand Internasional
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPEISIFICATION

Product Name

HEALGEN HIV 1/2 Human Immunodefficiency Virus Rapid Test

Reg. Number

KEMENKES RI AKL 30305917662

Specimen

Whole Blood / Serum / Plasma

Reading Time

15 Menit

Package

50 Test/Kit

Type

Cassette

Sensitivity

99,5 %

Spesificity

99,8 %

Function

The qualitative detection of antibodies to HIV-1 and/or HIV-2 in whole blood, serum or plasma

RAPID URINALYSIS TEST



URS-10 T (Urine Test Strip)

Reagent Strips for Urinalysis
Anti-VC Interference Ability



- ☒ Terdaftar di E-CATALOGUE LKPP
- ☒ Teregistrasi KEMENKES RI AKL
- ☒ Bersertifikasi Internasional
- ☒ Brand International
- ☒ Distributor bersertifikasi CDAKB



Link barcode produk:



SPESIFIKATION

Product Name

HEALGEN URS-10 T Reagent Strips for Urinalysis Anti-VC Interference Ability.

Reg. Number

KEMENKES RI AKL 20101718060

Specimen

Urine

Package

100 strips /Tube

Type

Strips

Parameter

Leukocytes, Nitrite, Urobilinogen, pH, Protein, Blood, Specific Gravity, Ketone, Bilirubin, Glucose



CFDA FDA



*e*catalogue



ISO : 9001

ISO : 13485