

WATO EX-55

The enhanced anesthesia workstation





Since 2006 Mindray has successfully installed over 10,000 units from the anesthesia range to customers all over the world. International feedback from our installed bases has reached an overwhelming level of customer satisfaction of which Mindray is immensely proud.

For the last decade, Mindray has continued to work closely with clinicians across the globe, to recognize and understand the clinical challenges encountered everyday and overcome them with new innovative and intuitive solutions. With this in mind Mindray is proud to bring you the enhanced WATO EX-55.

WATO EX-55

the enhanced anesthesia workstation

More Powerful

- Smarter ventilator technology
- Enhanced monitoring capabilities
- Complete range of ventilation mode

More Integrated

- Space saving design
- Plug-and-Play Multi-Gas module
- Built in optional patient suction system


More User-friendly

- 12.1 inch full color touch-screen
- Easy-to-use UI
- New flat-menu structure





More Powerful



Smarter Ventilator

The WATO EX-55 offers a smarter ventilator supporting a complete range of control and support ventilation modes to fully meet the needs of all your patients.



Driven Gas Auto-switch

To ensure the patient is ventilated at all times the WATO EX-55 has a unique driven gas auto-switch function. Should the primary ventilator driven gas be lost the device will automatically switch to Air (or O₂; dependent on device configuration) without interruption to mechanical ventilation.



Integrated, Heated Breathing System

With an integrated, user-controlled, heating system rain-out is reduced to a minimum providing improved flow sensor accuracy and guaranteeing extensive stability of tidal volume delivery. Both reusable and disposable soda lime canisters are compatible with the breathing circuit.





More Integrated

Compact Design, Smaller Footprint

The WATO EX-55 with a minimal footprint is ideal for the modern, busy and crowded OR suite. A space saving device that doesn't compromise on features or functionality.



PnP Multi-Gas Module

With improved single slot CO₂ module or double slot anesthetic agent capability the Mindray Plug-and-Play Multi-Gas modules provide comprehensive breath-by-breath analysis of FiO₂, EtO₂, CO₂, N₂O, auto-detection of five anesthetic agents.



Optional Integrated Suction for Better PeriOP Management

The new built in patient suction system (optional) provides the user with immediate access to a user-controlled suction system that is both neat and convenient.



More User-friendly



Intuitive UI Simplifies Your Workflow

Configured with a high resolution, intuitive, 12.1 inch full color touch-screen users are able to view and configure parameters as required. The intuitive layout and simple flat-menu structure ensure all parameters are clearly displayed and only two steps required for ventilation mode setting.



Large Easy-to-clean Worktable

The large, new, seamless and easy-to-clean writing table provides an optimized work surface for writing all your documents and keeping drugs and intubation equipment in a handy location.

