





BeneFusion DS3

Bedside Docking System

Data Sheet



Configuration

Description	Picture	Dimension	Remark
BeneFusion DS3 with 4 slots		719 x 162 x 79mm	Support 4 pumps; Use with Medical supply unit; Without pumps.
BeneFusion DS3 with 4 slots and four SP3		719 x 162 x 79mm	Use with Medical supply unit; With 4 SP3
BeneFusion DS3 with 6 slots, with four SP3 and two VP3, with trolley		1886 x 639 x 631mm	With four SP3 and two VP3;
BeneFusion DS3 with 6 slots, with four SP3 and two VP3		876 x 330 x 168mm	Use with Medical supply unit; With four SP3 and two VP3

Type

Classification Type CF, class I, IP22

Power supply

AC power 100-240V, 50/60Hz, 6-3A

Environment

Temperature 0-40°C for operating, -40-70°C for storage
Humidity 15-95% for operating, 10-95% for storage
Air pressure 57-106kPa for operating, 50-106kPa for storage

Mindray Building, Keji 12th Road South,
High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R.
China
Tel: +86 755 8188 8998 Fax: +86 755 26582680
E-mail: intl-market@mindray.com www.mindray.com

mindray | Healthcare with reach are registered trademarks or trademarks owned by Shenzhen Mindray Bio-medical Electronics Co., LTD
© 2020 Shenzhen Mindray Bio-medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.
P/N: ENG-BeneFusion DS3 Datasheet-210285x2P

mindray