

**PT MINDRAY MEDICAL INDONESIA**

**Mindray PM-60 Pulse Oximeter**

**Specifications**

Display	: 2.4 inch color TFT
Resolution:	: 320x240 pixels
Multi displays selectable	: Standard screen display, waveform display
Dimensions approx.	: 120 mmx55mmx30mm
Weight	: approx. max. 300g(with battery)
Trace	: 1 plethysmogram waveform
Battery	: 2 hours for charging, 24 hours for continuous working
Trend memory	: Resolution 2s, maximum time: 96 hours
	: Review of up to 4,000 data, measure up to 99 patients
Alarm	: User adjustable High and Low limits, 3- level audible and visual alarm

**SpO2**

Measurement range 0-100%

**Pulse Rate**

range 18 to 300bpm

Protection Against Ingress of Liquid

: IPX2

**Operation Environment**

Temperature	: 0 to 45 °C
Humidity	: 10% to 95%(non-condensing)
Temperature	: Operating: 0 to 45 °C
	: Storage temperature to – 20 °C)
Humidity	: Operating / Storage: 10 to 95%, non-condensation

**Included**

- : Main unit PM-60 with SpO2 & PR, 1pc
- : Charging kit (1 li-ion battery) & charger stand, power cable
- : Protective cover 1pc.