

MANUFACTURE STANDARD, MANUFACTURE RELIANCE

# RT-7200

## Auto Hematology Analyzer

NEW



### Features

- 20 parameters+3 histograms
- 3-part differentiation of WBC
- 2 Counting mode: whole blood and prediluted
- Throughput: 30 samples per hour
- 8.0" TFT color LCD with mouse operation
- Built-in thermal printer, optional external printer
- External AC adapter avoids electronic noises, enhancing accuracy
- Up to 30000 sample results (including histograms) can be stored
- Multi-language software available on request