



URIT-880

Chemistry Analyzer

- ✓ High-quality & close optical system
- ✓ Easy touch screen operation
- ✓ Convenient maintenance design
- Large memory capacity
- ✓ With 8 inner incubation positions
- ✓ Support barcode reader

URIT

URIT-880

Chemistry Analyzer

| Specifications |

Analysis method: Kinetic, end point, two point, dual wavelength, ABS

Light source: Halogen lamp, 6V/10W

Wavelength (nm): 340,405,492,510,546,578,630,700

Absorbance range: 0~4.0Abs

Wavelength accuracy: $\pm 2\text{nm}$

Spectrum band: $\pm 10\text{nm}$

Flowcell: 10mm quartz cell, 32 μL volume

Temperature control: Room temperature, 37°C ($\pm 0.1^{\circ}\text{C}$ variable)

Sampling volume: 100 μL ~ 9999 μL adjustable, recommended: $\geq 500 \mu\text{L}$

Cross contamination: $\leq 1.0\%$

Stability: $\leq 0.005\text{Abs/h}$

Memory: Up to 300 items and 30000 test results

Display: Large LCD touch screen display


Printer: Built-in thermal printer

Interface: Parallel port (For external printer), RJ-45, 3USB, SDID, PS-2


Ambient: Room temperature: 10°C ~ 32°C ; Humidity: $\leq 40\%$ ~ 85%

Power: AC (110-220)V, 50Hz

Dimension: 362.5 x 341 x 165mm (length * width * height)



7 inch touch screen



8 inner incubation positions

| Features |

- High-quality grating & close optical system
- Easy touch screen operation with 7 inch colour screen
- Convenient maintenance design
- Large memory capacity
- With 8 inner incubation positions
- Support barcode reader and keyboard

URIT MEDICAL ELECTRONIC CO., LTD.

Manufacturer: URIT Medical Electronic Co., Ltd.
Address: No. D-07 Information Industry District,
High-tech Zone, Guilin, Guangxi 541004,
P. R. China
Tel: +86-773-2288586 Fax: +86-773-2288560
E-mail: export@uritest.com service@uritest.com

Shenzhen Office
Add.: Room EJ, 24th Floor, New Baohui Building, Nanshan Road,
Nanshan District, Shenzhen, P.R. China
Tel: +86-755-26050853 Fax: +86-755-86036150
Http://www.urit.com

