



SARANA MAJU SEJAHTERA

VACUUM BLOOD COLLECTION TUBE



Item	Size (mm)	Volume	Tube Material	Additive	Color	Main Intended Use
Plain Tube	13x75	3 ml	Glass or PET	No Additive		Determinations in serum for clinical chemistry, Microbiological serology, Immunology
Pro-Coagulation Tube	13x75	3 ml	Glass or PET	Clot activator		
Gel & Clot Activator Tube	13x75 13x100	3 ml 4 ml 5 ml	Glass or PET	Gel & Clot Activator		
K3 EDTA Tube	10x45 13x75	0,5 ml 2 ml 3 ml	Glass or PET	K3 EDTA		Determinations in EDTA whole blood for hematology
PT Tube	13x75	1,8 ml 2,7 ml	PET	Sodium Citrate		Determinations in citrated plasma for coagulation testing
ESR Tube	8x120	1,28 ml	Glass	Sodium Citrate		Test for the blood cell sedimentation rate
Heparin Tube	13x75	3 ml	PET	Sodium Heparin or Lithium Heparin		Determinations in heparinized plasma for clinical chemistry
Glucose Tube	13x75	3 ml	PET	Inhibitor for Glucolysis Sodium Fluoride +EDTA Sodium Fluoride +Potassium oxalate		Determinations in stabilized anticoagulated whole blood or plasma for glucose and lactate testing